

Six1 Polyclonal Antibody

Cat No: HR1AP8317

For research use only

Overview

Product Name	Six1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	SIX1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SIX1; Homeobox protein SIX1; Sine oculis homeobox homolog 1
Background	SIX homeobox 1(SIX1) Homo sapiens The protein encoded by this gene is a homeobox protein that is similar to the Drosophila 'sine oculis' gene product. This gene is found in a cluster of related genes on chromosome 14 and is thought to be involved in limb development. Defects in this gene are a cause of autosomal dominant deafness type 23 (DFNA23) and branchiootic syndrome type 3 (BOS3). [provided by RefSeq, Jul 2008],