

Pbx 4 Polyclonal Antibody

Cat No: HR1AP8729

For research use only

Overview

Product Name	Pbx 4 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	PBX4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PBX4; Pre-B-cell leukemia transcription factor 4; Homeobox protein PBX4
Background	<p>PBX homeobox 4(PBX4) Homo sapiens This gene encodes a member of the pre-B cell leukemia transcription factor family. These proteins are homeobox proteins that play critical roles in embryonic development and cellular differentiation both as Hox cofactors and through Hox-independent pathways. The encoded protein contains a homeobox DNA-binding domain, but specific functions of the protein have not been determined. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, May 2011],</p>