

Nkx-6.1 Polyclonal Antibody

Cat No: HR1AP7505

For research use only

Overview

Product Name	Nkx-6.1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	NKX6-1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NKX6-1; NKX6A; Homeobox protein Nkx-6.1; Homeobox protein NK-6 homolog A
Background	NK6 homeobox 1(NKX6-1) Homo sapiens In the pancreas, NKX6.1 is required for the development of beta cells and is a potent bifunctional transcription regulator that binds to AT-rich sequences within the promoter region of target genes lype et al. (2004) [PubMed 15056733].[supplied by OMIM, Mar 2008],