

Nkx-2.4 Polyclonal Antibody

Cat No: HR1AP8289

For research use only

Overview

Product Name	Nkx-2.4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	NKX2-4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NKX2-4; NKX2D; Homeobox protein Nkx-2.4; Homeobox protein NK-2 homolog D
Background	function:Probable transcription factor.,similarity:Belongs to the NK-2 homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,