

Meis1 Polyclonal Antibody

Cat No: HR1AP7338

For research use only

Overview

Product Name	Meis1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	MEIS1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MEIS1; Homeobox protein Meis1
Background	<p>Meis homeobox 1(MEIS1) Homo sapiens Homeobox genes, of which the most well-characterized category is represented by the HOX genes, play a crucial role in normal development. In addition, several homeoproteins are involved in neoplasia. This gene encodes a homeobox protein belonging to the TALE (‘three amino acid loop extension’) family of homeodomain-containing proteins. [provided by RefSeq, Jul 2008],</p>