

# HoxA6 Polyclonal Antibody

Cat No: HR1AP3668

For research use only

## Overview

Product Name	HoxA6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	HOXA6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HOXA6; HOX1B; Homeobox protein Hox-A6; Homeobox protein Hox-1B
Background	homeobox A6(HOXA6) Homo sapiens In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation. [provided by RefSeq, Jul 2008],