

## Six2 Polyclonal Antibody

Cat No: HR1AP5658

For research use only

### Overview

Product Name	Six2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	SIX2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SIX2; Homeobox protein SIX2; Sine oculis homeobox homolog 2
Background	SIX homeobox 2(SIX2) Homo sapiens This gene is a member of the vertebrate gene family which encode proteins homologous to the Drosophila &apos;sine oculis&apos; homeobox protein. The encoded protein is a transcription factor which, like other members of this gene family, may be involved in limb or eye development. [provided by RefSeq, Dec 2008],