

Brn-3a Polyclonal Antibody

Cat No: HR1AP7774

For research use only

Overview

Product Name	Brn-3a Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	POU4F1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	POU4F1; BRN3A; RDC1; POU domain; class 4, transcription factor 1; Brain-specific homeobox/POU domain protein 3A; Brain-3A; Brn-3A; Homeobox/POU domain protein RDC-1; Oct-T1
Background	POU class 4 homeobox 1(POU4F1) Homo sapiens This gene encodes a member of the POU-IV class of neural transcription factors. This protein is expressed in a subset of retinal ganglion cells and may be involved in the developing sensory nervous system. This protein may also promote the growth of cervical tumors. A translocation of this gene is associated with some adult acute myeloid leukemias. [provided by RefSeq, Mar 2012],