

## Nkx-6.3 Polyclonal Antibody

Cat No: HR1AP6056

For research use only

### Overview

Product Name	Nkx-6.3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	NKX6-3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NKX6-3; Homeobox protein Nkx-6.3
Background	NK6 homeobox 3(NKX6-3) Homo sapiens The NKX family of homeodomain proteins controls numerous developmental processes. Members of the NKX6 subfamily, including NKX6-3, are involved in development of the central nervous system (CNS), gastrointestinal tract, and pancreas (Alanentalo et al., 2006 [PubMed 16326147]).[supplied by OMIM, Mar 2008],