

Dlx-5 Polyclonal Antibody

Cat No: HR1AP3289

For research use only

Overview

Product Name	Dlx-5 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GeneID?Human?	DLX5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DLX5; Homeobox protein DLX-5
Background	<p>distal-less homeobox 5(DLX5) Homo sapiens This gene encodes a member of a homeobox transcription factor gene family similiar to the Drosophila distal-less gene. The encoded protein may play a role in bone development and fracture healing. Mutation in this gene, which is located in a tail-to-tail configuration with another member of the family on the long arm of chromosome 7, may be associated with split-hand/split-foot malformation. [provided by RefSeq, Jul 2008],</p>