

# ALX3 Polyclonal Antibody

Cat No: HR1AP2755

For research use only

## Overview

Product Name	ALX3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	ALX3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ALX3; Homeobox protein aristaless-like 3; Proline-rich transcription factor ALX3
Background	ALX homeobox 3(ALX3) Homo sapiens This gene encodes a nuclear protein with a homeobox DNA-binding domain that functions as a transcriptional regulator involved in cell-type differentiation and development. Preferential methylation of this gene's promoter is associated with advanced-stage neuroblastoma tumors. [provided by RefSeq, Jul 2008],