

# PREP Polyclonal Antibody

Cat No: HR1AP4973

For research use only

## Overview

Product Name	PREP Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	PREP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PREP; PEP; Prolyl endopeptidase; PE; Post-proline cleaving enzyme
Background	<p>prolyl endopeptidase(PREP) Homo sapiens The protein encoded by this gene is a cytosolic prolyl endopeptidase that cleaves peptide bonds on the C-terminal side of prolyl residues within peptides that are up to approximately 30 amino acids long. Prolyl endopeptidases have been reported to be involved in the maturation and degradation of peptide hormones and neuropeptides. [provided by RefSeq, Jul 2008],</p>