

GHITM Polyclonal Antibody

Cat No: HR1AP6648

For research use only

Overview

Product Name	GHITM Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	GHITM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GHITM; DERP2; MICS1; TMBIM5; My021; Growth hormone-inducible transmembrane protein; Dermal papilla-derived protein 2; Mitochondrial morphology and cristae structure 1; MICS1; Transmembrane BAX inhibit
Background	similarity:Belongs to the BI1 family.,