

CTRP8 Polyclonal Antibody

Cat No: HR1AP7184

For research use only

Overview

Product Name	CTRP8 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	C1QTNF8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	C1QTNF8; Complement C1q tumor necrosis factor-related protein 8; C1q/TNF-related protein 8; CTRP8
Background	similarity:Contains 1 C1q domain.,similarity:Contains 1 collagen-like domain.,