

Inhibin β -C Polyclonal Antibody

Cat No: HR1AP9516

For research use only

Overview

Product Name	Inhibin β -C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	INHBC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	INHBC; Inhibin beta C chain; Activin beta-C chain
Background	<p>inhibin beta C subunit(INHBC) Homo sapiens This gene encodes a member of the TGF-beta (transforming growth factor-beta) superfamily of proteins. The encoded preproprotein is proteolytically processed to generate a subunit of homodimeric and heterodimeric activin complexes. The heterodimeric complex may function in the inhibition of activin A signaling. Transgenic mice overexpressing this gene exhibit defects in testis, liver and prostate. [provided by RefSeq, Aug 2016].</p>