

IL-26 Polyclonal Antibody

Cat No: HR1AP9513

For research use only

Overview

Product Name	IL-26 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	IL26
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IL26; AK155; Interleukin-26; IL-26; Protein AK155
Background	interleukin 26(IL26) Homo sapiens This gene was identified by its overexpression specifically in herpesvirus samimiri-transformed T cells. The encoded protein is a member of the IL10 family of cytokines. It is a secreted protein and may function as a homodimer. This protein is thought to contribute to the transformed phenotype of T cells after infection by herpesvirus samimiri. [provided by RefSeq, Jul 2008],