

IFN- ϵ Polyclonal Antibody

Cat No: HR1AP5145

For research use only

Overview

Product Name	IFN- ϵ Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	IFNE
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IFNE; IFNE1; Interferon epsilon; IFN-epsilon; Interferon epsilon-1
Background	similarity:Belongs to the alpha/beta interferon family.,tissue specificity:Expressed in prostate cancer cell line PC-3, amnion-derived WISH cells, SK-MEL28 melanoma cells and Daudi cells and very weakly in MCF-7 human breast cancer cells. Not detected in Namalwa cells (Burkitt's lymphoma) or in Jurkat (T) cells.,