

# INSL4 Polyclonal Antibody

Cat No: HR1AP3735

For research use only

## Overview

Product Name	INSL4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GeneID?Human?	INSL4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	INSL4; Early placenta insulin-like peptide; EPIL; Insulin-like peptide 4; Placentin
Background	insulin like 4(INSL4) Homo sapiens INSL4 encodes the insulin-like 4 protein, a member of the insulin superfamily. INSL4 encodes a precursor that undergoes post-translational cleavage to produce 3 polypeptide chains, A-C, that form tertiary structures composed of either all three chains, or just the A and B chains. Expression of INSL4 products occurs within the early placental cytotrophoblast and syncytiotrophoblast. [provided by RefSeq, Jul 2008],