

## IL-20 Polyclonal Antibody

Cat No: HR1AP7576

For research use only

### Overview

Product Name	IL-20 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	IL20
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IL20; ZCYTO10; Interleukin-20; IL-20; Cytokine Zcyto10
Background	interleukin 20(IL20) Homo sapiens The protein encoded by this gene is a cytokine structurally related to interleukin 10 (IL10). This cytokine has been shown to transduce its signal through signal transducer and activator of transcription 3 (STAT3) in keratinocytes. A specific receptor for this cytokine is found to be expressed in skin and upregulated dramatically in psoriatic skin, suggesting a role for this protein in epidermal function and psoriasis. [provided by RefSeq, Jul 2008],