

Cardiotrophin-1 Polyclonal Antibody

Cat No: HR1AP2953

For research use only

Overview

Product Name	Cardiotrophin-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	CTF1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CTF1; Cardiotrophin-1; CT-1
Background	cardiotrophin 1(CTF1) Homo sapiens The protein encoded by this gene is a secreted cytokine that induces cardiac myocyte hypertrophy in vitro. It has been shown to bind and activate the ILST/gp130 receptor. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2008],