

IL-20 Polyclonal Antibody

Cat No: HR1AP9830

For research use only

Overview

Product Name	IL-20 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	IL20 ZCYTO10 UNQ852/PRO1801
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Interleukin-20 (IL-20) (Cytokine Zcyto10)
Background	interleukin 20(IL20) Homo sapiens The protein encoded by this gene is a cytokine structurally related to interleukin 10 (IL10). This cytokine has been shown to transduce its signal through signal transducer and activator of transcription 3 (STAT3) in keratinocytes. A specific receptor for this cytokine is found to be expressed in skin and upregulated dramatically in psoriatic skin, suggesting a role for this protein in epidermal function and psoriasis. [provided by RefSeq, Jul 2008].