

## Epigen Polyclonal Antibody

Cat No: HR1AP5492

For research use only

### Overview

Product Name	Epigen Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GeneID?Human?	EPGN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EPGN; Epigen; Epithelial mitogen; EPG
Background	epithelial mitogen(EPGN) Homo sapiens The protein encoded by this gene is a member of the epidermal growth factor family. Members of this family are ligands for the epidermal growth factor receptor and play a role in cell survival, proliferation and migration. This protein has been reported to have high mitogenic activity but low affinity for its receptor. Expression of this transcript and protein have been reported in cancer specimens of the breast, bladder, and prostate. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012],