

FGF-16 Polyclonal Antibody

Cat No: HR1AP5282

For research use only

Overview

Product Name	FGF-16 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	FGF16
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FGF16; Fibroblast growth factor 16; FGF-16
Background	fibroblast growth factor 16(FGF16) Homo sapiens This gene encodes a member of a family of proteins that are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This gene is expressed in cardiac cells and is required for proper heart development. Mutation in this gene was also observed in individuals with metacarpal 4-5 fusion. [provided by RefSeq, Mar 2014],