

GMF-? Polyclonal Antibody

Cat No: HR1AP5354

For research use only

Overview

Product Name	GMF-? Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GeneID?Human?	GMFB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GMFB; Glia maturation factor beta; GMF-beta
Background	function:This protein causes differentiation of brain cells, stimulation of neural regeneration, and inhibition of proliferation of tumor cells.,PTM:Phosphorylated; stimulated by phorbol ester.,similarity:Belongs to the actin-binding proteins ADF family. GMF subfamily.,similarity:Contains 1 ADF-H domain.,