

IL-22 Polyclonal Antibody

Cat No: HR1AP9512

For research use only

Overview

Product Name	IL-22 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	IL22
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IL22; ILTIF; ZCYTO18; Interleukin-22; IL-22; Cytokine Zcyto18; IL-10-related T-cell-derived-inducible factor; IL-TIF
Background	function:Cytokine that contributes to the inflammatory response in vivo.,similarity:Belongs to the IL-10 family.,