

Cripto-3 Polyclonal Antibody

Cat No: HR1AP9511

For research use only

Overview

Product Name	Cripto-3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	TDGF1P3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TDGF1P3; CRIPTO3; TDGF2; TDGF3; Putative teratocarcinoma-derived growth factor 3; Cripto-3 growth factor; Epidermal growth factor-like cripto protein CR3; Teratocarcinoma-derived growth factor 1 pseud
Background	teratocarcinoma-derived growth factor 1 pseudogene 3(TDGF1P3) Homo sapiens The TDGF3 locus has characteristics of a retrotransposon, including lack of introns and a poly(A) sequence. [provided by RefSeq, Jul 2008],