

IFN- ω Polyclonal Antibody

Cat No: HR1AP5256

For research use only

Overview

Product Name	IFN- ω Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	IFNW1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IFNW1; Interferon omega-1; Interferon alpha-II-1
Background	interferon omega 1(IFNW1) Homo sapiens The protein encoded by this gene is an interferon and possesses antiviral activity. The encoded protein binds to the interferon alpha/beta receptor but not to the interferon gamma receptor. This intronless gene has several pseudogenes spread throughout the genome. [provided by RefSeq, Nov 2015],