

HRP-3 Polyclonal Antibody

Cat No: HR1AP5325

For research use only

Overview

Product Name	HRP-3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	HDGFRP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HDGFRP3; HDGF2; CGI-142; Hepatoma-derived growth factor-related protein 3; HRP-3; Hepatoma-derived growth factor 2; HDGF-2
Background	function:Enhances DNA synthesis and may play a role in cell proliferation.,similarity:Belongs to the HDGF family.,similarity:Contains 1 PWWP domain.,tissue specificity:Detected in testis, heart, spinal cord and brain.,