

# Gastrin Polyclonal Antibody

Cat No: HR1AP6521

For research use only

## Overview

Product Name	Gastrin Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	GAST
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GAST; GAS; Gastrin
Background	gastrin(GAST) Homo sapiens Gastrin is a hormone whose main function is to stimulate secretion of hydrochloric acid by the gastric mucosa, which results in gastrin formation inhibition. This hormone also acts as a mitogenic factor for gastrointestinal epithelial cells. Gastrin has two biologically active peptide forms, G34 and G17. [provided by RefSeq, Jul 2008],