

ET-1 Polyclonal Antibody

Cat No: HR1AP4988

For research use only

Overview

Product Name	ET-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	EDN1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EDN1; Endothelin-1; Preoendothelin-1; PPET1
Background	<p>endothelin 1(EDN1) Homo sapiens This gene encodes a preproprotein that is proteolytically processed to generate a secreted peptide that belongs to the endothelin/sarafotoxin family. This peptide is a potent vasoconstrictor and its cognate receptors are therapeutic targets in the treatment of pulmonary arterial hypertension. Aberrant expression of this gene may promote tumorigenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2015],</p>