

HGFA Polyclonal Antibody

Cat No: HR1AP4911

For research use only

Overview

Product Name	HGFA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	HGFAC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HGFAC; Hepatocyte growth factor activator; HGF activator; HGFA
Background	HGF activator(HGFAC) Homo sapiens This gene encodes a member of the peptidase S1 protein family. The encoded protein is first synthesized as an inactive single-chain precursor before being activated to a heterodimeric form by endoproteolytic processing. It acts as serine protease that converts hepatocyte growth factor to the active form. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014],