

Prolactin Polyclonal Antibody

Cat No: HR1AP5365

For research use only

Overview

Product Name	Prolactin Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	PRL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PRL; Prolactin; PRL
Background	<p>prolactin(PRL) Homo sapiens This gene encodes the anterior pituitary hormone prolactin. This secreted hormone is a growth regulator for many tissues, including cells of the immune system. It may also play a role in cell survival by suppressing apoptosis, and it is essential for lactation. Alternative splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Aug 2011],</p>