

ADM Polyclonal Antibody

Cat No: HR1AP5043

For research use only

Overview

Product Name	ADM Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	ADM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ADM; AM; ADM
Background	<p>adrenomedullin(ADM) Homo sapiens The protein encoded by this gene is a preprohormone which is cleaved to form two biologically active peptides, adrenomedullin and proadrenomedullin N-terminal 20 peptide. Adrenomedullin is a 52 aa peptide with several functions, including vasodilation, regulation of hormone secretion, promotion of angiogenesis, and antimicrobial activity. The antimicrobial activity is antibacterial, as the peptide has been shown to kill <i>E. coli</i> and <i>S. aureus</i> at low concentration. [provided by RefSeq, Aug 2014].</p>