

Sg III Polyclonal Antibody

Cat No: HR1AP4983

For research use only

Overview

Product Name	Sg III Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	53kDa
GeneID?Human?	SCG3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SCG3; Secretogranin-3; Secretogranin III; SgIII
Background	secretogranin III(SCG3) Homo sapiens The protein encoded by this gene is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. Granins may serve as precursors for biologically active peptides. Some granins have been shown to function as helper proteins in sorting and proteolytic processing of prohormones; however, the function of this protein is unknown. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009],