

Inhibin β -B Polyclonal Antibody

Cat No: HR1AP5249

For research use only

Overview

Product Name	Inhibin β -B Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	INHBB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	INHBB; Inhibin beta B chain; Activin beta-B chain
Background	<p>inhibin beta B subunit(INHBB) Homo sapiens This gene encodes a member of the TGF-beta (transforming growth factor-beta) superfamily of proteins. The encoded preproprotein is proteolytically processed to generate a subunit of the dimeric activin and inhibin protein complexes. These complexes activate and inhibit, respectively, follicle stimulating hormone secretion from the pituitary gland. Polymorphisms near this gene are associated with pre-eclampsia in female human patients. [provided by RefSeq, Aug 2016],</p>