

Shd Polyclonal Antibody

Cat No: HR1AP7975

For research use only

Overview

Product Name	Shd Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	130kDa
GeneID?Human?	SHD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SHD; SH2 domain-containing adapter protein D
Background	function:May function as an adapter protein.,PTM:Tyrosine phosphorylated by ABL.,similarity:Contains 1 SH2 domain.,