

# TAGAP Polyclonal Antibody

Cat No: HR1AP5781

For research use only

## Overview

Product Name	TAGAP Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	TAGAP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TAGAP; TAGAP1; FKSG15; T-cell activation Rho GTPase-activating protein; T-cell activation GTPase-activating protein
Background	T-cell activation RhoGTPase activating protein(TAGAP) Homo sapiens This gene encodes a member of the Rho GTPase-activator protein superfamily. The encoded protein may function as a Rho GTPase-activating protein. Alterations in this gene may be associated with several diseases, including rheumatoid arthritis, celiac disease, and multiple sclerosis. Alternate splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2013],