

## RERG Polyclonal Antibody

Cat No: HR1AP8906

For research use only

### Overview

Product Name	RERG Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	RERG
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RERG; Ras-related and estrogen-regulated growth inhibitor
Background	RAS like estrogen regulated growth inhibitor(RERG) Homo sapiens RERG, a member of the RAS superfamily of GTPases, inhibits cell proliferation and tumor formation (Finlin et al., 2001 [PubMed 11533059]).[supplied by OMIM, Mar 2009],