

# RASSF7 Polyclonal Antibody

Cat No: HR1AP8728

For research use only

## Overview

Product Name	RASSF7 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	RASSF7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RASSF7; C11orf13; HRC1; Ras association domain-containing protein 7; HRAS1-related cluster protein 1
Background	similarity:Contains 1 Ras-associating domain.,