

# RASSF2 Polyclonal Antibody

Cat No: HR1AP9209

For research use only

## Overview

Product Name	RASSF2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	RASSF2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RASSF2; KIAA0168; Ras association domain-containing protein 2
Background	Ras association domain family member 2(RASSF2) Homo sapiens This gene encodes a protein that contains a Ras association domain. Similar to its cattle and sheep counterparts, this gene is located near the prion gene. Two alternatively spliced transcripts encoding the same isoform have been reported. [provided by RefSeq, Jul 2008],