

RASL10B Polyclonal Antibody

Cat No: HR1AP9086

For research use only

Overview

Product Name	RASL10B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	RASL10B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RASL10B; Ras-like protein family member 10B; Ras-like protein VTS58635; Ras-related protein 17; RRP17
Background	similarity:Belongs to the small GTPase superfamily. Ras family.,