

## Rap 2B Polyclonal Antibody

Cat No: HR1AP8101

For research use only

### Overview

Product Name	Rap 2B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	RAP2B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAP2B; Ras-related protein Rap-2b
Background	<p>RAP2B, member of RAS oncogene family(RAP2B) Homo sapiens This intronless gene belongs to a family of RAS-related genes. The proteins encoded by these genes share approximately 50% amino acid identity with the classical RAS proteins and have numerous structural features in common. The most striking difference between the RAP and RAS proteins resides in their 61st amino acid: glutamine in RAS is replaced by threonine in RAP proteins. Evidence suggests that this protein may be polyisoprenylated and palmitoylated. [provided by RefSeq, Jul 2008],</p>