

Rab 5C Polyclonal Antibody

Cat No: HR1AP8076

For research use only

Overview

Product Name	Rab 5C Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	RAB5C
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAB5C; RABL; Ras-related protein Rab-5C; L1880; RAB5L
Background	RAB5C, member RAS oncogene family(RAB5C) Homo sapiens Members of the Rab protein family are small GTPases of the Ras superfamily that are thought to ensure fidelity in the process of docking and/or fusion of vesicles with their correct acceptor compartment (Han et al., 1996 [PubMed 8646882]).[supplied by OMIM, Nov 2010],