

## Rab 37 Polyclonal Antibody

Cat No: HR1AP6885

For research use only

### Overview

Product Name	Rab 37 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	RAB37
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAB37; Ras-related protein Rab-37
Background	RAB37, member RAS oncogene family(RAB37) Homo sapiens Rab proteins are low molecular mass GTPases that are critical regulators of vesicle trafficking. For additional background information on Rab proteins, see MIM 179508.[supplied by OMIM, Apr 2006],