

Rab 2B Polyclonal Antibody

Cat No: HR1AP8899

For research use only

Overview

Product Name	Rab 2B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	RAB2B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAB2B; Ras-related protein Rab-2B
Background	RAB2B, member RAS oncogene family(RAB2B) Homo sapiens Members of the Rab protein family are nontransforming monomeric GTP-binding proteins of the Ras superfamily that contain 4 highly conserved regions involved in GTP binding and hydrolysis. Rab proteins are prenylated, membrane-bound proteins involved in vesicular fusion and trafficking; see MIM 179508.[supplied by OMIM, Apr 2006],