

# Rab 24 Polyclonal Antibody

Cat No: HR1AP7741

For research use only

## Overview

Product Name	Rab 24 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	RAB24
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAB24; Ras-related protein Rab-24
Background	RAB24, member RAS oncogene family(RAB24) Homo sapiens RAB24 is a small GTPase of the Rab subfamily of Ras-related proteins that regulate intracellular protein trafficking (Olikkonen et al., 1993 [PubMed 8126105]).[supplied by OMIM, Aug 2009],