

Rab 11B Polyclonal Antibody

Cat No: HR1AP9107

For research use only

Overview

Product Name	Rab 11B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	RAB11B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RAB11B; YPT3; Ras-related protein Rab-11B; GTP-binding protein YPT3
Background	RAB11B, member RAS oncogene family(RAB11B) Homo sapiens The Ras superfamily of small GTP-binding proteins, which includes the Ras (see MIM 190020), Ral (see MIM 179550), Rho (see MIM 165390), Rap (see MIM 179520), and Rab (see MIM 179508) families, is involved in controlling a diverse set of essential cellular functions. The Rab family, including RAB11B, appears to play a critical role in regulating exocytotic and endocytotic pathways (summary by Zhu et al., 1994 [PubMed 7811277]).[supplied by OMIM, Nov 2010],