

## GTPBP5 Polyclonal Antibody

Cat No: HR1AP6580

For research use only

### Overview

Product Name	GTPBP5 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	GTPBP5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GTPBP5; OBGH1; GTP-binding protein 5; Protein obg homolog 1; ObgH1
Background	mitochondrial ribosome associated GTPase 2(MTG2) Homo sapiens Small G proteins, such as GTPBP5, act as molecular switches that play crucial roles in the regulation of fundamental cellular processes such as protein synthesis, nuclear transport, membrane trafficking, and signal transduction (Hirano et al., 2006 [PubMed 17054726]).[supplied by OMIM, Mar 2008],