

GTPBP2 Polyclonal Antibody

Cat No: HR1AP7778

For research use only

Overview

Product Name	GTPBP2 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	GTPBP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GTPBP2; GTP-binding protein 2
Background	GTP binding protein 2(GTPBP2) Homo sapiens GTP-binding proteins, or G proteins, constitute a superfamily capable of binding GTP or GDP. G proteins are activated by binding GTP and are inactivated by hydrolyzing GTP to GDP. This general mechanism enables G proteins to perform a wide range of biologic activities.[supplied by OMIM, Jan 2003],