

EF-Tu Polyclonal Antibody

Cat No: HR1AP8573

For research use only

Overview

Product Name	EF-Tu Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	TUFM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TUFM; Elongation factor Tu; mitochondrial; EF-Tu; P43
Background	Tu translation elongation factor, mitochondrial(TUFM) Homo sapiens This gene encodes a protein which participates in protein translation in mitochondria. Mutations in this gene have been associated with combined oxidative phosphorylation deficiency resulting in lactic acidosis and fatal encephalopathy. A pseudogene has been identified on chromosome 17. [provided by RefSeq, Jul 2008],